

PTBP2 Blocking Peptide for STJ506038 peptide (STJ505790)

STJ505790

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PTBP2 Blocking Peptide for STJ506038 is synthetically produced from the 150-200 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	58155
Gene Symbol	PTBP2
Uniprot ID	PTBP2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 150-200 on human Polypyrimidine tract-binding protein 2/PTBP2.
Immunogen Region	150-200
Specificity	This blocking peptide is recommended for use in combination with PTBP2 antibody, STJ506038
Immunogen Sequence	