

TUFM Positive Control for STJ506048 peptide (STJ505797)

STJ505797

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	TUFM Positive Control for STJ506048 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7284
Gene Symbol	TUFM
Uniprot ID	EFTU_HUMAN
Immunogen	Purified protein.
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with TUFM antibody STJ506048.
Immunogen Sequence	