

PI3R6 Positive Control for STJ506052 peptide (STJ505799)

STJ505799

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PI3R6 Positive Control for STJ506052 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	146850
Gene Symbol	PIK3R6
Uniprot ID	PI3R6_HUMAN
Immunogen	Purified protein.
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with PI3R6 antibody STJ506052.
Immunogen Sequence	