

PP1R3A Blocking Peptide for STJ506059 peptide (STJ505804)

STJ505804

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	PP1R3A Blocking Peptide for STJ506059 is synthetically produced from the 830-880 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5506
Gene Symbol	PPP1R3A
Uniprot ID	PPP1R3A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 830-880 on mouse Protein Phosphatase 1 Regulatory Subunit 3A.
Immunogen Region	830-880
Specificity	This blocking peptide is recommended for use in combination with PP1R3A antibody, STJ506059
Immunogen Sequence	