

## ARHGEF3 Blocking Peptide peptide (Phospho) (STJ505836)

STJ505836

### GENERAL INFORMATION

<b>Product Type</b>	Biomolecules
<b>Short Description</b>	Phospho-ARHGEF3 Blocking Peptide is synthetically produced from the 350-410 sequence and is suitable for use in western blot applications.
<b>Applications</b>	Immunodepletion/Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	50650
<b>Gene Symbol</b>	ARHGEF3
<b>Uniprot ID</b>	ARHG3_HUMAN
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 350-410 human Rho-Guanine nucleotide Exchange factor 3 protein.
<b>Immunogen Region</b>	350-410
<b>Specificity</b>	
<b>Immunogen Sequence</b>	