

ARHGEF3 Positive Control peptide (Phospho) (STJ505837)

STJ505837

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-ARHGEF3 Positive Control is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50650
Gene Symbol	ARHGEF3
Uniprot ID	ARHG3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	