

CGRP Blocking Peptide for STJ506112 peptide (STJ505838)

STJ505838

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	CGRP Blocking Peptide for STJ506112 is synthetically produced from the 350-410 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	12310
Gene Symbol	Calca
Uniprot ID	CALCA_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 350-410 mouse Calcitonin gene-related peptide 1.
Immunogen Region	350-410
Specificity	This blocking peptide is recommended for use in combination with CGRP antibody, STJ506112
Immunogen Sequence	