

Cofilin Blocking Peptide peptide (Phospho) (STJ505844)

STJ505844

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-Cofilin Blocking Peptide is synthetically produced from the 1-50 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1072
Gene Symbol	CFL1
Uniprot ID	COF1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 1-50 on human Phosphorylated Cofilin.
Immunogen Region	1-50
Specificity	
Immunogen	
Sequence	