

C1ORF115 Blocking Peptide for STJ506143 peptide (STJ505863)

STJ505863

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	C1ORF115 Blocking Peptide for STJ506143 is synthetically produced from the 50-100 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	226777
Gene Symbol	
Uniprot ID	RDD1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 50-100 on mouse Chromosome 1 Open Reading Frame 115.
Immunogen Region	50-100
Specificity	This blocking peptide is recommended for use in combination with C1ORF115 antibody, STJ506143
Immunogen Sequence	