

C1ORF115 Positive Control for STJ506143 peptide (STJ505880)

STJ505880

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	C1ORF115 Positive Control for STJ506143 is synthetically produced from the 75-125 sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	79762
Gene Symbol	C1orf115
Uniprot ID	RDD1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 75-125 on human Chromosome 1 Open Reading Frame 115.
Immunogen Region	75-125
Specificity	This is positive control is recommended for use in combination with C1ORF115 antibody STJ506143.
Immunogen Sequence	