

NEDD8 Positive Control for STJ506209 peptide (STJ505901)

STJ505901

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	NEDD8 Positive Control for STJ506209 is synthetically produced from the sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4738
Gene Symbol	NEDD8
Uniprot ID	NEDD8_HUMAN
Immunogen	
Immunogen Region	
Specificity	This is positive control is recommended for use in combination with NEDD8 antibody STJ506209.
Immunogen Sequence	