

Anti-PP2A-beta antibody (450-497) (STJ505947)

STJ505947

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PP2A-beta (450-497) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IHC/IP/WB
Host/Source	Rabbit
Reactivity	Human/Monkey/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5526
Gene Symbol	PPP2R5B
Uniprot ID	2A5B_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-497 on human Serine/Threonine-Protein Phosphatase 2A 56 kDa regulatory Beta subunit.
Immunogen Region	450-497
Specificity	
Immunogen Sequence	