

**Anti-PP2A-beta antibody (450-497) {FITC} (STJ505949)**

STJ505949

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-PP2A-beta (450-497) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA/IP/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Monkey/Mouse/Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:200
<b>Formulation</b>	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">5526</a>
<b>Gene Symbol</b>	<a href="#">PPP2R5B</a>
<b>Uniprot ID</b>	<a href="#">2A5B_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 450-497 on human Serine/Threonine-Protein Phosphatase 2A 56 kDa regulatory Beta subunit.
<b>Immunogen Region</b>	450-497
<b>Specificity</b>	
<b>Immunogen Sequence</b>	