

Anti-PP2A-delta antibody (540-600) {Biotin} (STJ505951)

STJ505951

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PP2A-delta (540-600) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5528
Gene Symbol	PPP2R5D
Uniprot ID	2A5D_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 540-600 on human Serine/Threonine-Protein Phosphatase 2A B56-Delta subunit.
Immunogen Region	540-600
Specificity	
Immunogen Sequence	