

Anti-CXCL5 antibody (1-60) (STJ506004)

STJ506004

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-CXCL5 (1-60) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA/IP/WB |
| Host/Source | Rabbit |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.5-1.5 µg/µl |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 20311 |
| Gene Symbol | Cxcl5 |
| Uniprot ID | CXCL5_MOUSE |
| Immunogen | Synthetic peptide taken within amino acid region 1-60 on mouse C-X-C motif chemokine 5. |
| Immunogen Region | 1-60 |
| Specificity | |
| Immunogen Sequence | |