

Anti-Phospho-PBXIP1 antibody (100-150) {Biotin} (STJ506021)

STJ506021

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-PBXIP1 (100-150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	229534
Gene Symbol	Pbxip1
Uniprot ID	PBIP1_MOUSE
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 100-150 on mouse Pre-B-cell leukemia transcription factor-interacting protein 1.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	