

Anti-PI3R6 antibody (590-640) {FITC} (STJ506052)

STJ506052

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PI3R6 (590-640) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146850
Gene Symbol	PIK3R6
Uniprot ID	PI3R6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 590-640 on human Phosphoinositide 3-kinase regulatory subunit 6.
Immunogen Region	590-640
Specificity	
Immunogen Sequence	