

Anti-ASXL2 antibody (1200-1275) (STJ506147)

STJ506147

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ASXL2 (1200-1275) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	75302
Gene Symbol	Asxl2
Uniprot ID	ASXL2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1200-1275 on mouse Putative Polycomb group protein ASXL2.
Immunogen Region	1200-1275
Specificity	
Immunogen Sequence	