

Anti-Phospho-PDE11A-S163 antibody (100-200) {Biotin} (STJ506174)

STJ506174

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-PDE11A-S163 (100-200) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	140928
Gene Symbol	Pde11a
Uniprot ID	PDE11_RAT
Immunogen	Synthetic peptide taken within amino acid region 100-200 on rat PDE11A protein. Phosphorylated Ser 163.
Immunogen Region	100-200
Specificity	
Immunogen Sequence	