

Anti-RNase H1 antibody (240-286) {FITC} (STJ506226)

STJ506226

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RNase H1 (240-286) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	246243
Gene Symbol	RNASEH1
Uniprot ID	RNH1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 240-286 on human Rnase H1 protein.
Immunogen Region	240-286
Specificity	
Immunogen Sequence	