

**Anti-SLC26A6 antibody (700-759) (STJ506236)**

STJ506236

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-SLC26A6 (700-759) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA/IP/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 Mu g/Mu l
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100 ICC: 1:50-1:100 IF: 1:50-1:100
<b>Formulation</b>	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">65010</a>
<b>Gene Symbol</b>	<a href="#">SLC26A6</a>
<b>Uniprot ID</b>	<a href="#">S26A6_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 700-759 on human SCL26A6 protein.
<b>Immunogen Region</b>	700-759
<b>Specificity</b>	
<b>Immunogen Sequence</b>	