

Dnmt1 Blocking Peptide for STJ506245 peptide (STJ506246)

STJ506246

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Dnmt1 Blocking Peptide for STJ506245 is synthetically produced from the 1550-1600 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1786
Gene Symbol	DNMT1
Uniprot ID	DNMT1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 1550-1600 on human Dnmt1 protein.
Immunogen Region	1550-1600
Specificity	This blocking peptide is recommended for use in combination with Dnmt1 antibody, STJ506245
Immunogen Sequence	