

Dnmt2 Blocking Peptide for STJ506249 peptide (STJ506250)

STJ506250

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Dnmt2 Blocking Peptide for STJ506249 is synthetically produced from the 350-400 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	13434
Gene Symbol	Trdmt1
Uniprot ID	TRDMT_MOUSE
Immunogen	350-400
Immunogen Region	350-400
Specificity	This blocking peptide is recommended for use in combination with Dnmt2 antibody, STJ506249
Immunogen Sequence	