

Anti-Dnmt3a antibody (100-150) (STJ506252)

STJ506252

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dnmt3a (100-150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	ELISA 1:10, 000 ICC 1:50-1:100 IF 1:50-1:100 IHC 1:50-1:100 IP 1:200 WB 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13435
Gene Symbol	Dnmt3a
Uniprot ID	DNM3A_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 100-150 on mouse DNA (cytosine-5)-methyltransferase 3A protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	