

Anti-Dnmt3L antibody (375-421) {FITC} (STJ506264)

STJ506264

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dnmt3L (375-421) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	ELISA 1:10, 000 ICC 1:50-1:100 IF 1:50-1:100 IHC 1:50-1:100 IP 1:200 WB 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54427
Gene Symbol	Dnmt3l
Uniprot ID	DNM3L_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 375-421 on mouse DNA (cytosine-5)-methyltransferase 3-like protein.
Immunogen Region	375-421
Specificity	
Immunogen Sequence	