

NEDD4 Blocking Peptide for STJ506278 peptide (STJ506279)

STJ506279

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	NEDD4 Blocking Peptide for STJ506278 is synthetically produced from the 300-350 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4734
Gene Symbol	NEDD4
Uniprot ID	NEDD4_HUMAN
Immunogen	300-350
Immunogen Region	300-350
Specificity	This blocking peptide is recommended for use in combination with NEDD4 antibody, STJ506278
Immunogen Sequence	