

**Anti-ATP citrate synthase Catalytic Subunit antibody (600-650) {Biotin} (STJ506309)**

STJ506309

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-ATP citrate synthase Catalytic Subunit (600-650) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA/IP/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Yeast

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 Mu g/Mu l
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	ELISA 1:10, 000 ICC 1:50-1:100 IF 1:50-1:100 IHC 1:50-1:100 IP 1:200 WB 1:500
<b>Formulation</b>	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 600-650 on yeast ATP Citrate Lyase Catalytic subunit protein.
<b>Immunogen Region</b>	600-650
<b>Specificity</b>	
<b>Immunogen Sequence</b>	