

Amyloid Precursor Protein Positive Control for STJ506314 peptide (STJ506318)

STJ506318

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Amyloid Precursor Protein Positive Control for STJ506314 is synthetically produced from the 700-770 sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	351
Gene Symbol	APP
Uniprot ID	A4_HUMAN
Immunogen	700-770
Immunogen Region	700-770
Specificity	This is positive control is recommended for use in combination with Amyloid Precursor Protein antibody STJ506314.
Immunogen Sequence	