

Anti-CEACAM1 antibody (480-526) {FITC} (STJ506321)

STJ506321

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CEACAM1 (480-526) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	ELISA 1:10, 000 IP 1:200 WB 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	634
Gene Symbol	CEACAM1
Uniprot ID	CEAM1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 480-526 on human Carcinoembryonic antigen-related Cell Adhesion Molecule 1 protein.
Immunogen Region	480-526
Specificity	
Immunogen Sequence	