

Anti-Skn1 antibody (150-200) {FITC} (STJ506336)

STJ506336

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Skn1 (150-200) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	ELISA 1:10, 000 IP 1:200 WB 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25833
Gene Symbol	POU2F3
Uniprot ID	P02F3_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 150-200 on human POU2F3 protein.
Immunogen Region	150-200
Specificity	
Immunogen Sequence	