

Keap1 Positive Control for STJ506354 peptide (STJ506358)

STJ506358

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Biomolecules |
| Short Description | Keap1 Positive Control for STJ506354 is synthetically produced from the 500-550 sequence and is suitable for use in western blot applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 9817 |
| Gene Symbol | KEAP1 |
| Uniprot ID | KEAP1_HUMAN |
| Immunogen | 500-550 |
| Immunogen Region | 500-550 |
| Specificity | This is positive control is recommended for use in combination with Keap1 antibody STJ506354. |
| Immunogen Sequence | |