

TREM2 Positive Control for STJ506374 peptide (STJ506378)

STJ506378

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	TREM2 Positive Control for STJ506374 is synthetically produced from the 100-200 sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9365
Gene Symbol	KL
Uniprot ID	KL0T_HUMAN
Immunogen	100-200
Immunogen Region	100-200
Specificity	This is positive control is recommended for use in combination with TREM2 antibody STJ506374.
Immunogen Sequence	