

SLCO6C1 Positive Control for STJ506398 peptide (STJ506402)

STJ506402

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	SLCO6C1 Positive Control for STJ506398 is synthetically produced from the 500-600 sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	74441
Gene Symbol	Slco6c1
Uniprot ID	SO6C1_MOUSE
Immunogen	500-600
Immunogen Region	500-600
Specificity	This is positive control is recommended for use in combination with SLCO6C1 antibody STJ506398.
Immunogen Sequence	