

Anti-Phospho-SPT5 antibody (760-810) {FITC} (STJ506405)

STJ506405

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-SPT5 (760-810) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	ELISA 1:10, 000 IP 1:200 WB 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20924
Gene Symbol	Supt5h
Uniprot ID	SPT5H_MOUSE
Immunogen	Phosphorylated synthetic peptide corresponding to amino acids within region 760-810 on mouse SPT5 protein.
Immunogen Region	760-810
Specificity	
Immunogen Sequence	