

SPT5 Blocking Peptide for STJ506404 peptide (Phospho) (STJ506406)

STJ506406

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Phospho-SPT5 Blocking Peptide for STJ506404 is synthetically produced from the 760-810 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	20924
Gene Symbol	Supt5h
Uniprot ID	SPT5H_MOUSE
Immunogen	760-810
Immunogen Region	760-810
Specificity	This blocking peptide is recommended for use in combination with SPT5 antibody, STJ506404
Immunogen Sequence	