

ADCY2 Blocking Peptide for STJ506425 peptide (STJ506426)

STJ506426

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	ADCY2 Blocking Peptide for STJ506425 is synthetically produced from the C-Term sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	IP 1:200
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	Adcy2
Uniprot ID	ADCY2_MOUSE
Immunogen	C-epitope
Immunogen Region	C-Term
Specificity	This blocking peptide is recommended for use in combination with ADCY2 antibody, STJ506425
Immunogen Sequence	