

Anti-GCDH/GCD antibody (200-250) {FITC} (STJ507758)

STJ507758

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GCDH/GCD (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2639
Gene Symbol	GCDH
Uniprot ID	GCDH_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 200-250 on Human GCDH/GCD protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	