

Anti-SKOR1 antibody (C-Term) {FITC} (STJ507763)

STJ507763

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SKOR1 (C-Term) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24473
Gene Symbol	Htr1a
Uniprot ID	5HT1A_RAT
Immunogen	Synthetic peptide corresponding to amino acids within region 850-900 on Human SKOR1 protein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	