

SKOR1 Positive Control for STJ507768 peptide (STJ507770)

STJ507770

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	SKOR1 Positive Control for STJ507768 is synthetically produced from the N-Term sequence and is suitable for use in western blot applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	24473
Gene Symbol	Htr1a
Uniprot ID	5HT1A_RAT
Immunogen	Synthetic peptide corresponding to amino acids within region 650-700 on human SKOR1 protein.
Immunogen Region	N-Term
Specificity	This is positive control is recommended for use in combination with SKOR1 antibody STJ507768.
Immunogen Sequence	