

Anti-ADCY2 antibody (Multiple Epitope) {Biotin} (STJ507775)

STJ507775

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-ADCY2 (Multiple Epitope) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	Adcy2
Uniprot ID	ADCY2_MOUSE
Immunogen	Synthetic peptide corresponding to peptides taken from amino acid region 850-900 and 950-1000 on mouse ADCY2 protein.
Immunogen Region	Multiple Epitope
Specificity	
Immunogen Sequence	