

Anti-Lysosomal acid lipase antibody (200-250) {Biotin} (STJ507785)

STJ507785

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lysosomal acid lipase (200-250) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3988
Gene Symbol	LIPA
Uniprot ID	LICH_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 200-250 on Human Lysosomal acid lipase protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	