

Adeno-Associated Virus Capsid Protein Blocking Peptide for STJ507812 peptide (STJ507815)

STJ507815

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	Adeno-Associated Virus Capsid Protein Blocking Peptide for STJ507812 is synthetically produced from the Mid Region sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	Virus

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to amino acids within region 200-250 on rat 5-hydroxytryptamine receptor 1A protein.
Immunogen Region	Mid Region
Specificity	This blocking peptide is recommended for use in combination with Adeno-Associated Virus Capsid Protein antibody, STJ507812
Immunogen Sequence	