

Anti-SATB2 antibody (400-450) {FITC} (STJ507822)

STJ507822

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SATB2 (400-450) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23314
Gene Symbol	SATB2
Uniprot ID	SATB2_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 400-450 on human DNA-binding protein SATB2 protein.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	