

Anti-Dopamine beta hydroxylase antibody (125-178) {Biotin} (STJ507831)

STJ507831

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dopamine beta hydroxylase (125-178) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	13166
Gene Symbol	Dbh
Uniprot ID	DOPQ_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 125-178 on mouse Dopamine beta hydroxylase.
Immunogen Region	125-178
Specificity	
Immunogen Sequence	