

Anti-PLAP antibody (C-Term) {Biotin} (STJ507853)

STJ507853

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Placental alkaline phosphatase (PLAP) (C-Term) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA/IP/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	250
Gene Symbol	ALPP
Uniprot ID	PPB1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 480-535 on human Placental alkaline phosphatase (PLAP).
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	