

SALL4 Blocking Peptide for STJ507860 peptide (STJ507862)

STJ507862

GENERAL INFORMATION

Product Type	Biomolecules
Short Description	SALL4 Blocking Peptide for STJ507860 is synthetically produced from the 400-450 sequence and is suitable for use in western blot applications.
Applications	Immunodepletion/Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57167
Gene Symbol	SALL4
Uniprot ID	SALL4_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 400-450 on human DNA-binding protein SATB2 protein.
Immunogen Region	400-450
Specificity	This blocking peptide is recommended for use in combination with SALL4 antibody, STJ507860
Immunogen Sequence	