

Anti-MTMR1-Isoform 1 antibody (C-Term) (STJ70084)

STJ70084

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-MTMR1-Isoform 1 (C-Term) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave no signal but low background in Human Liver and Placenta lysates at up to 1µg/ml. We would appreciate any feedback from people in the field-have any result
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8776
Gene Symbol	MTMR1
Uniprot ID	MTMR1_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	Please note this will only recognise Human isoform 1 (NP_003819.1) , not isoform 2 (NP_789746).
Immunogen Sequence	SSPSHSATSVHTSV