

Anti-Dysadherin antibody (C-Term) (STJ70141)

STJ70141

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Dysadherin (C-Term) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA/WB
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. WB: Approx 38kDa band observed in Human Spleen lysates (calculated MW of 19.5 kDa according to NP_054883 and NP_659003). The band was successfully blocked by incubation with the immunizing peptide.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53827
Gene Symbol	FXVD5
Uniprot ID	FXVD5_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	NP_054883.3, NP_659003.1 and NP_001158077.1 are variants that represent the same protein.
Immunogen Sequence	GKCRQLSRLCRNHCR



STJ70141 staining (0.5 µg/ml) of Human Spleen lysate (RIPA buffer, 30 µg total protein per lane). Primary incubated for 1 hour. Detected by chemiluminescence.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081