

Anti-FOXD2 antibody (C-Term) (STJ70387)

STJ70387

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-FOXD2 (C-Term) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. WB: No signal obtained yet but low background observed in Human Kidney and 293 lysates at up to 1µg/ml. We would appreciate any feedback from people in the field-have any results been reported
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2306
Gene Symbol	FOXD2
Uniprot ID	FOXD2_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	FASKVAGLSGCHF