

Anti-RANBP16/Exportin 7 antibody (C-Term) (STJ70494)

STJ70494

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-RANBP16/Exportin 7 (C-Term) is suitable for use in ELISA, Western Blot and Immunohistochemistry research applications.
Applications	Pep-ELISA/WB/IHC
Host/Source	Goat
Reactivity	Human/Mouse/Dog/Pig/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 110kDa band observed in Hela lysates (predicted size of approx. 124kDa according to NP_055839). Recommended concentration: 0.5-1 µg/ml IHC: Paraffin embedded Human Testis. Re
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage	Store at -20°C on receipt and minimise freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	23039
Gene Symbol	XPO7
Uniprot ID	XPO7_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	This antibody is expected to recognise all the reported isoforms (NP_001093631.1; NP_055839; NP_001093632.1).
Immunogen Sequence	NSTYGVNSNDMMMS



STJ70494 staining (0.5 µg/ml) of Hela lysate (RIPA buffer, 35 µg total protein per lane). Primary incubated for 1 hour. Detected by chemiluminescence.